

Antibodies Product Datasheet

Product Name:	HSP105 (Phospho-Ser809) Antibody
Catalog Number:	A9092
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human HSP105 (Phospho-Ser809)
Concentration:	1 mg/ml
Gene Name:	HSPH1 HSP105 HSP110 KIAA0201
Protein Name:	HSP105 (Phospho-Ser809)
Other Name:	Heat shock protein 105 kDa (Antigen NY-CO-25) (Heat shock 110 kDa protein)
Human Swiss Prot No.:	Q92598
Formulation:	Liquid in PBS containing 50% glycerol, and 0.218% sodium azide.
Specificity:	This antibody detects endogenous levels of HSP105 (Phospho-Ser809) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year