

Antibodies Product Datasheet

Product Name:	Myd88 (Phospho-Tyr257) Antibody
Catalog Number:	A9099
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Myd88 (Phospho-Tyr257)
Concentration:	1 mg/ml
Gene Name:	MYD88
Protein Name:	Myd88 (Phospho-Tyr257)
Other Name:	Myeloid differentiation primary response protein MyD88
Human Swiss Prot No.:	Q99836
Formulation:	Liquid in PBS containing 50% glycerol, and 0.102% sodium azide.
Specificity:	This antibody detects endogenous levels of Myd88 (Phospho-Tyr257) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year