

Antibodies Product Datasheet

Product Name:	NDEL1 (Phospho-Thr219) Antibody
Catalog Number:	A9101
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human NDEL1 (Phospho-Thr219)
Concentration:	1 mg/ml
Gene Name:	NDEL1 EOPA MITAP1 NUDEL
Protein Name:	NDEL1 (Phospho-Thr219)
Other Name:	Nuclear distribution protein nudE-like 1 (Protein Nudel) (Mitosis-associated protein 1)
Human Swiss Prot No.:	Q9GZM8
Formulation:	Liquid in PBS containing 50% glycerol, and 0.219% sodium azide.
Specificity:	This antibody detects endogenous levels of NDEL1 (Phospho-Thr219) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year