

Antibodies Product Datasheet

Product Name:	Catenin δ -1 (Phospho-Tyr228) Antibody
Catalog Number:	A9402
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human Catenin δ -1 (Tyr228)
Concentration:	1 mg/ml
Gene Name:	CTNND1 KIAA0384
Protein Name:	Catenin δ -1 (Tyr228)
Other Name:	Catenin delta-1 (Cadherin-associated Src substrate) (CAS) (p120 catenin) (p120(ctn)) (p120(cas))
Human Swiss Prot No.:	O60716
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Catenin δ -1 (phospho-Tyr228)
Storage and Stability:	-20°C/1 year