

Antibodies Product Datasheet

Product Name:	cdc2 (Phospho-Thr14) Antibody
Catalog Number:	A9407
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human cdc2 (Thr14)
Concentration:	1 mg/ml
Gene Name:	CDK1 CDC2 CDC28A CDKN1 P34CDC2
Protein Name:	cdc2 (Thr14)
Other Name:	Cyclin-dependent kinase 1 (CDK1) (EC 2.7.11.22) (EC 2.7.11.23) (Cell division control protein 2 homolog) (Cell division protein kinase 1) (p34 protein kinase)
Human Swiss Prot No.:	P06493
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human cdc2 (phospho-Thr14)
Storage and Stability:	-20°C/1 year