

Antibodies Product Datasheet

Product Name:	CDC20 (Phospho-Ser51) Antibody
Catalog Number:	A9408
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human CDC20 (Ser51)
Concentration:	1 mg/ml
Gene Name:	CDC20
Protein Name:	CDC20 (Ser51)
Other Name:	Cell division cycle protein 20 homolog (p55CDC)
Human Swiss Prot No.:	Q12834
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse CDC20 (phospho-Ser51)
Storage and Stability:	-20°C/1 year