

Antibodies Product Datasheet

Product Name:	Chk2 (Phospho-Thr432) Antibody
Catalog Number:	A9411
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Chk2 (Phospho Thr432)
Concentration:	1 mg/ml
Gene Name:	CHEK2 CDS1 CHK2 RAD53
Protein Name:	Chk2 (Phospho Thr432)
Other Name:	Serine/threonine-protein kinase Chk2 (EC 2.7.11.1) (CHK2 checkpoint homolog) (Cds1 homolog) (Hucds1) (hCds1) (Checkpoint kinase 2)
Human Swiss Prot No.:	O96017
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat Chk2 (Phospho Thr432)
Storage and Stability:	-20°C/1 year