

Antibodies Product Datasheet

Product Name:	Cyclin B1 (Phospho-Ser116) Antibody
Catalog Number:	A9413
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human Cyclin B1 (Ser116)
Concentration:	1 mg/ml
Gene Name:	CCNB1 CCNB
Protein Name:	Cyclin B1 (Ser116)
Other Name:	G2/mitotic-specific cyclin-B1
Human Swiss Prot No.:	P14635
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Cyclin B1 (phospho-Ser116)
Storage and Stability:	-20°C/1 year