

Antibodies Product Datasheet

Product Name:	Cyclin D1 (Phospho-Ser90) Antibody
Catalog Number:	A9415
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Cyclin D1 (Phospho Ser90)
Concentration:	1 mg/ml
Gene Name:	CCND1 BCL1 PRAD1
Protein Name:	Cyclin D1 (Phospho Ser90)
Other Name:	G1/S-specific cyclin-D1 (B-cell lymphoma 1 protein) (BCL-1) (BCL-1 oncogene) (PRAD1 oncogene)
Human Swiss Prot No.:	P24385
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat Cyclin D1 (Phospho Ser90)
Storage and Stability:	-20°C/1 year