

Antibodies Product Datasheet

Product Name:	Cyclin E1 (Phospho-Thr62) Antibody
Catalog Number:	A9416
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human Cyclin E1 (Thr62)
Concentration:	1 mg/ml
Gene Name:	CCNE1 CCNE
Protein Name:	Cyclin E1 (Thr62)
Other Name:	G1/S-specific cyclin-E1
Human Swiss Prot No.:	P24864
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Cyclin E1 (phospho-Thr62)
Storage and Stability:	-20°C/1 year