

Antibodies Product Datasheet

Product Name:	GABBR2 (Phospho-Ser893) Antibody
Catalog Number:	A9423
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human GABBR2 (Phospho Ser893)
Concentration:	1 mg/ml
Gene Name:	GABBR2 GPR51 GPRC3B
Protein Name:	GABBR2 (Phospho Ser893)
Other Name:	Gamma-aminobutyric acid type B receptor subunit 2 (GABA-B receptor 2) (GABA-B-R2) (GABA-BR2) (GABABR2) (Gb2) (G-protein coupled receptor 51) (HG20)
Human Swiss Prot No.:	O75899
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat GABBR2 (Phospho Ser893)
Storage and Stability:	-20°C/1 year