

Antibodies Product Datasheet

Product Name:	IP3 Receptor (Phospho-Tyr353) Antibody
Catalog Number:	A9424
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human IP3 Receptor (Phospho Tyr353)
Concentration:	1 mg/ml
Gene Name:	ITPR1 INSP3R1
Protein Name:	IP3 Receptor (Phospho Tyr353)
Other Name:	Inositol 1,4,5-trisphosphate receptor type 1 (IP3 receptor isoform 1) (IP3R 1) (InsP3R1) (Type 1 inositol 1,4,5-trisphosphate receptor) (Type 1 InsP3 receptor)
Human Swiss Prot No.:	Q14643
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat IP3 Receptor (Phospho Tyr353)
Storage and Stability:	-20°C/1 year