

Antibodies Product Datasheet

Product Name:	LIMK-2 (Phospho-Thr505) Antibody
Catalog Number:	A9429
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized phospho-peptide around the phosphorylation site of human LIMK-2 (phospho Thr505)
Concentration:	1 mg/ml
Gene Name:	LIMK2
Protein Name:	LIM domain kinase 2
Other Name:	LIMK2; LIM domain kinase 2; LIMK-2
Human Swiss Prot No.:	P53671
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-LIMK-2 (T505) Polyclonal Antibody detects endogenous levels of LIMK-2 protein only when phosphorylated at T505.
Storage and Stability:	-20°C/1 year