

Antibodies Product Datasheet

Product Name:	ATM (Phospho-Ser1987) Antibody
Catalog Number:	A9434
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ATM (Phospho Ser1987)
Concentration:	1 mg/ml
Gene Name:	ATM
Protein Name:	ATM (Phospho Ser1987)
Other Name:	Serine-protein kinase ATM (EC 2.7.11.1) (Ataxia telangiectasia mutated) (A-T mutated)
Human Swiss Prot No.:	Q62388(Q13315)
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat ATM (Phospho Ser1987)
Storage and Stability:	-20°C/1 year