

Antibodies Product Datasheet

Product Name:	ATR (Phospho-Thr1989) Antibody
Catalog Number:	A9435
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human ATR (Thr1989)
Concentration:	1 mg/ml
Gene Name:	ATR FRP1
Protein Name:	ATR (Thr1989)
Other Name:	Serine/threonine-protein kinase ATR (EC 2.7.11.1) (Ataxia telangiectasia and Rad3-related protein) (FRAP-related protein 1)
Human Swiss Prot No.:	Q13535
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human ATR (phospho-Thr1989)
Storage and Stability:	-20°C/1 year