

Antibodies Product Datasheet

Product Name:	Axl (Phospho-Tyr702) Antibody
Catalog Number:	A9436
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human Axl (Tyr702)
Concentration:	1 mg/ml
Gene Name:	AURKA AIK AIRK1 ARK1 AURA AYK1 BTAK IAK1 STK15 STK6
Protein Name:	Axl (Tyr702)
Other Name:	Aurora kinase A (EC 2.7.11.1) (Aurora 2) (Aurora/IPL1-related kinase 1) (ARK-1) (Aurora-related kinase 1) (hARK1) (Breast tumor-amplified kinase) (Serine/threonine-protein kinase 15) (Serine/threonine-protein kinase 6) (Serine/threonine-protein kinase aur
Human Swiss Prot No.:	O14965
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Axl (phospho-Tyr702)
Storage and Stability:	-20°C/1 year