

Antibodies Product Datasheet

Product Name:	Btk (Phospho-Ser180) Antibody
Catalog Number:	A9437
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human Btk (Ser180)
Concentration:	1 mg/ml
Gene Name:	BTK AGMX1 ATK BPK
Protein Name:	Btk (Ser180)
Other Name:	Tyrosine-protein kinase BTK (EC 2.7.10.2) (Agammaglobulinaemia tyrosine kinase) (ATK) (B-cell progenitor kinase) (BPK) (Bruton tyrosine kinase)
Human Swiss Prot No.:	Q06187
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Btk (phospho-Ser180)
Storage and Stability:	-20°C/1 year