

Antibodies Product Datasheet

Product Name:	c-Abl (Phospho-Tyr89) Antibody
Catalog Number:	A9441
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human c-Abl (Tyr89)
Concentration:	1 mg/ml
Gene Name:	ABL1 ABL JTK7
Protein Name:	c-Abl (Tyr89)
Other Name:	Tyrosine-protein kinase ABL1 (EC 2.7.10.2) (Abelson murine leukemia viral oncogene homolog 1) (Abelson tyrosine-protein kinase 1) (Proto-oncogene c-Abl) (p150)
Human Swiss Prot No.:	P00519
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human c-Abl (phospho-Tyr89)
Storage and Stability:	-20°C/1 year