

Antibodies Product Datasheet

Product Name:	CaMKK2 (Phospho-Ser511) Antibody
Catalog Number:	A9442
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human CaMKK2 (Ser511)
Concentration:	1 mg/ml
Gene Name:	CAMKK2 CAMKKB KIAA0787
Protein Name:	CaMKK2 (Ser511)
Other Name:	Calcium/calmodulin-dependent protein kinase kinase 2 (CaM-KK 2) (CaM-kinase kinase 2) (CaMKK 2) (EC 2.7.11.17) (Calcium/calmodulin-dependent protein kinase kinase beta) (CaM-KK beta) (CaM-kinase kinase beta) (CaMKK beta)
Human Swiss Prot No.:	Q96RR4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Rat CaMKK2 (phospho-Ser511)
Storage and Stability:	-20°C/1 year