

Antibodies Product Datasheet

Product Name:	GSK-3 β (Phospho-Thr390) Antibody
Catalog Number:	A9444
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human GSK-3 β (Thr390)
Concentration:	1 mg/ml
Gene Name:	GSK3B
Protein Name:	GSK-3 β (Thr390)
Other Name:	Glycogen synthase kinase-3 beta (GSK-3 beta) (EC 2.7.11.26) (Serine/threonine-protein kinase GSK3B) (EC 2.7.11.1)
Human Swiss Prot No.:	P49841
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human GSK-3 β (phospho-Thr390)
Storage and Stability:	-20°C/1 year