

Antibodies Product Datasheet

Product Name:	IKK ϵ (Phospho-Ser172) Antibody
Catalog Number:	A9445
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human IKK ϵ (Ser172)
Concentration:	1 mg/ml
Gene Name:	IKBKE IKKE IKKI KIAA0151
Protein Name:	IKK ϵ (Ser172)
Other Name:	Inhibitor of nuclear factor kappa-B kinase subunit epsilon (I-kappa-B kinase epsilon) (IKK-E) (IKK-epsilon) (IkbKE) (EC 2.7.11.10) (Inducible I kappa-B kinase) (IKK-i)
Human Swiss Prot No.:	Q14164
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human IKK ϵ (phospho-Ser172)
Storage and Stability:	-20°C/1 year