

Antibodies Product Datasheet

Product Name:	Jak1(Phospho-Tyr1034/1035) Antibody
Catalog Number:	A9446
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human Jak1(Tyr1034 and 1035)
Concentration:	1 mg/ml
Gene Name:	JAK1 JAK1A JAK1B
Protein Name:	Jak1(Tyr1034/1035)
Other Name:	Tyrosine-protein kinase JAK1 (EC 2.7.10.2) (Janus kinase 1) (JAK-1)
Human Swiss Prot No.:	P23458
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Rat Jak1(phospho-Tyr1034 or 1035)
Storage and Stability:	-20°C/1 year