

Antibodies Product Datasheet

Product Name:	Jak2 (Phospho-Tyr1008) Antibody
Catalog Number:	A9447
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human Jak2 (Tyr1008)
Concentration:	1 mg/ml
Gene Name:	JAK2
Protein Name:	Jak2 (Tyr1008)
Other Name:	Tyrosine-protein kinase JAK2 (EC 2.7.10.2) (Janus kinase 2) (JAK-2)
Human Swiss Prot No.:	O60674
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Jak2 (phospho-Tyr1008)
Storage and Stability:	-20°C/1 year