

Antibodies Product Datasheet

Product Name:	MKK7 (Phospho-Ser271/Thr275) Antibody
Catalog Number:	A9452
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human MKK7 (Ser271 and Thr275)
Concentration:	1 mg/ml
Gene Name:	MAP2K7 JNKK2 MEK7 MKK7 PRKMK7 SKK4
Protein Name:	MKK7 (Ser271/Thr275)
Other Name:	Dual specificity mitogen-activated protein kinase kinase 7 (MAP kinase kinase 7) (MAPKK 7) (EC 2.7.12.2) (JNK-activating kinase 2) (MAPK/ERK kinase 7) (MEK 7) (Stress-activated protein kinase kinase 4) (SAPK kinase 4) (SAPKK-4) (SAPKK4) (c-Jun N-terminal
Human Swiss Prot No.:	O14733
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human MKK7 (phospho-Ser271 or Thr275)
Storage and Stability:	-20°C/1 year