

Antibodies Product Datasheet

Product Name:	CD79A (Phospho-Tyr182) Antibody
Catalog Number:	A9487
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human CD79A (Tyr182)
Concentration:	1 mg/ml
Gene Name:	CD79A IGA MB1
Protein Name:	CD79A (Tyr182)
Other Name:	B-cell antigen receptor complex-associated protein alpha chain (Ig-alpha) (MB-1 membrane glycoprotein) (Membrane-bound immunoglobulin-associated protein) (Surface IgM-associated protein) (CD antigen CD79a)
Human Swiss Prot No.:	P11912
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human CD79A (phospho-Tyr182)
Storage and Stability:	-20°C/1 year