

Antibodies Product Datasheet

Product Name:	LAT (Phospho-Tyr171) Antibody
Catalog Number:	A9490
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human LAT (Tyr171)
Concentration:	1 mg/ml
Gene Name:	LAT
Protein Name:	LAT (Tyr171)
Other Name:	Linker for activation of T-cells family member 1 (36 kDa phospho-tyrosine adapter protein) (pp36) (p36-38)
Human Swiss Prot No.:	O43561
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human LAT (phospho-Tyr171)
Storage and Stability:	-20°C/1 year