

## Antibodies Product Datasheet

|                               |                                                                                                                    |
|-------------------------------|--------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | EGF Receptor (Phospho-Tyr992) Antibody                                                                             |
| <b>Catalog Number:</b>        | A9502                                                                                                              |
| <b>Reactivity:</b>            | Human, Mouse                                                                                                       |
| <b>Applications:</b>          | WB                                                                                                                 |
| <b>Dilutions:</b>             | WB 1:1000-2000                                                                                                     |
| <b>Modification:</b>          | Phospho Specific                                                                                                   |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                  |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.          |
| <b>Immunogen:</b>             | Synthesized phospho peptide around human EGF Receptor (Tyr992)                                                     |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                            |
| <b>Gene Name:</b>             | EGFR ERBB ERBB1 HER1                                                                                               |
| <b>Protein Name:</b>          | EGF Receptor (Tyr992)                                                                                              |
| <b>Other Name:</b>            | Epidermal growth factor receptor (EC 2.7.10.1) (Proto-oncogene c-ErbB-1) (Receptor tyrosine-protein kinase erbB-1) |
| <b>Human Swiss Prot No.:</b>  | P00533                                                                                                             |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                            |
| <b>Specificity:</b>           | This antibody detects endogenous levels of Human Mouse EGF Receptor (phospho-Tyr992)                               |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                       |