

Antibodies Product Datasheet

Product Name:	ALK (Phospho-Tyr1278/1282/1283) Antibody
Catalog Number:	A9505
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human ALK (Tyr1278 and 1282 and 1283)
Concentration:	1 mg/ml
Gene Name:	ALK
Protein Name:	ALK (Tyr1278/1282/1283)
Other Name:	ALK tyrosine kinase receptor (EC 2.7.10.1) (Anaplastic lymphoma kinase) (CD antigen CD246)
Human Swiss Prot No.:	Q9UM73
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human ALK (phospho-Tyr1278 or 1282 or 1283)
Storage and Stability:	-20°C/1 year