

Antibodies Product Datasheet

Product Name:	c-Kit (Phospho-Tyr568/570) Antibody
Catalog Number:	A9508
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human c-Kit (Tyr568 and 570)
Concentration:	1 mg/ml
Gene Name:	KIT SCFR
Protein Name:	c-Kit (Tyr568/570)
Other Name:	Mast/stem cell growth factor receptor Kit (SCFR) (EC 2.7.10.1) (Piebald trait protein) (PBT) (Proto-oncogene c-Kit) (Tyrosine-protein kinase Kit) (p145 c-kit) (v-kit Hardy-Zuckerman 4 feline sarcoma viral oncogene homolog) (CD antigen CD117)
Human Swiss Prot No.:	P10721
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human c-Kit (phospho-Tyr568 or 570)
Storage and Stability:	-20°C/1 year