

Antibodies Product Datasheet

Product Name:	DDR1 (Phospho-Tyr796) Antibody
Catalog Number:	A9511
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DDR1 (Phospho Tyr796)
Concentration:	1 mg/ml
Gene Name:	DDR1 CAK EDDR1 NEP NTRK4 PTK3A RTK6 TRKE
Protein Name:	DDR1 (Phospho Tyr796)
Other Name:	Epithelial discoidin domain-containing receptor 1 (Epithelial discoidin domain receptor 1) (EC 2.7.10.1) (CD167 antigen-like family member A) (Cell adhesion kinase) (Discoidin receptor tyrosine kinase) (HGK2) (Mammary carcinoma kinase 10) (MCK-10) (Protein)
Human Swiss Prot No.:	Q08345
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat DDR1 (Phospho Tyr796)
Storage and Stability:	-20°C/1 year