

Antibodies Product Datasheet

Product Name:	DDR1 (Phospho-Tyr792) Antibody
Catalog Number:	A9512
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human DDR1 (Tyr792)
Concentration:	1 mg/ml
Gene Name:	DDR1 CAK EDDR1 NEP NTRK4 PTK3A RTK6 TRKE
Protein Name:	DDR1 (Tyr792)
Other Name:	Epithelial discoidin domain-containing receptor 1 (Epithelial discoidin domain receptor 1) (EC 2.7.10.1) (CD167 antigen-like family member A) (Cell adhesion kinase) (Discoidin receptor tyrosine kinase) (HGK2) (Mammary carcinoma kinase 10) (MCK-10) (Protein)
Human Swiss Prot No.:	Q08345
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human DDR1 (phospho-Tyr792)
Storage and Stability:	-20°C/1 year