

Antibodies Product Datasheet

Product Name:	EphA2 (Phospho-Tyr772) Antibody
Catalog Number:	A9515
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human EphA2 (Tyr772)
Concentration:	1 mg/ml
Gene Name:	EPHA2 ECK
Protein Name:	EphA2 (Tyr772)
Other Name:	Ephrin type-A receptor 2 (EC 2.7.10.1) (Epithelial cell kinase) (Tyrosine-protein kinase receptor ECK)
Human Swiss Prot No.:	P29317
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human EphA2 (phospho-Tyr772)
Storage and Stability:	-20°C/1 year