

Antibodies Product Datasheet

Product Name:	EPHA3 (Phospho-Tyr602) Antibody
Catalog Number:	A9516
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EPHA3 (Phospho Tyr602)
Concentration:	1 mg/ml
Gene Name:	EPHA3 ETK ETK1 HEK TYRO4
Protein Name:	EPHA3 (Phospho Tyr602)
Other Name:	Ephrin type-A receptor 3 (EC 2.7.10.1) (EPH-like kinase 4) (EK4) (hEK4) (HEK) (Human embryo kinase) (Tyrosine-protein kinase TYRO4) (Tyrosine-protein kinase receptor ETK1) (Eph-like tyrosine kinase 1)
Human Swiss Prot No.:	P29320
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat EPHA3 (Phospho Tyr602)
Storage and Stability:	-20°C/1 year