

Antibodies Product Datasheet

Product Name:	FGF Receptor (Phospho-Tyr653/654) Antibody
Catalog Number:	A9519
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human FGF Receptor (Tyr653 and 654)
Concentration:	1 mg/ml
Gene Name:	FGFR1 BFGFR CEK FGFBR FLG FLT2 HBGFR
Protein Name:	FGF Receptor (Tyr653/654)
Other Name:	Fibroblast growth factor receptor 1 (FGFR-1) (EC 2.7.10.1) (Basic fibroblast growth factor receptor 1) (BFGFR) (bFGF-R-1) (Fms-like tyrosine kinase 2) (FLT-2) (N-sam) (Proto-oncogene c-Fgr) (CD antigen CD331)
Human Swiss Prot No.:	P11362
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human FGF Receptor (phospho-Tyr653 or 654)
Storage and Stability:	-20°C/1 year