

Antibodies Product Datasheet

Product Name:	IGF-I Receptor β (Phospho-Tyr980) Antibody
Catalog Number:	A9528
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human IGF-I Receptor β (Tyr980)
Concentration:	1 mg/ml
Gene Name:	IGF1R
Protein Name:	IGF-I Receptor β (Tyr980)
Other Name:	Insulin-like growth factor 1 receptor (EC 2.7.10.1) (Insulin-like growth factor I receptor) (IGF-I receptor) (CD antigen CD221) [Cleaved into: Insulin-like growth factor 1 receptor alpha chain; Insulin-like growth factor 1 receptor beta chain]
Human Swiss Prot No.:	P08069
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Rat IGF-I Receptor β (phospho-Tyr980)
Storage and Stability:	-20°C/1 year