

Antibodies Product Datasheet

Product Name:	Insulin Receptor β (Phospho-Tyr1345) Antibody
Catalog Number:	A9529
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human Insulin Receptor β (Tyr1345)
Concentration:	1 mg/ml
Gene Name:	INSR
Protein Name:	Insulin Receptor β (Tyr1345)
Other Name:	Insulin receptor (IR) (EC 2.7.10.1) (CD antigen CD220) [Cleaved into: Insulin receptor subunit alpha; Insulin receptor subunit beta]
Human Swiss Prot No.:	P06213
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Insulin Receptor β (phospho-Tyr1345)
Storage and Stability:	-20°C/1 year