

Antibodies Product Datasheet

Product Name:	M-CSF Receptor (Phospho-Tyr546) Antibody
Catalog Number:	A9531
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human M-CSF Receptor (Tyr546)
Concentration:	1 mg/ml
Gene Name:	CSF1R FMS
Protein Name:	M-CSF Receptor (Tyr546)
Other Name:	Macrophage colony-stimulating factor 1 receptor (CSF-1 receptor) (CSF-1-R) (CSF-1R) (M-CSF-R) (EC 2.7.10.1) (Proto-oncogene c-Fms) (CD antigen CD115)
Human Swiss Prot No.:	P07333
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse M-CSF Receptor (phospho-Tyr546)
Storage and Stability:	-20°C/1 year