

Antibodies Product Datasheet

Product Name:	Met (Phospho-Tyr1234/1235) Antibody
Catalog Number:	A9534
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human Met (Tyr1234 and 1235)
Concentration:	1 mg/ml
Gene Name:	MET
Protein Name:	Met (Tyr1234/1235)
Other Name:	Hepatocyte growth factor receptor (HGF receptor) (EC 2.7.10.1) (HGF/SF receptor) (Proto-oncogene c-Met) (Scatter factor receptor) (SF receptor) (Tyrosine-protein kinase Met)
Human Swiss Prot No.:	P08581
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Rat Met (phospho-Tyr1234 or 1235)
Storage and Stability:	-20°C/1 year