

Antibodies Product Datasheet

Product Name:	APP (Phospho-Thr668) Antibody
Catalog Number:	A9543
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human APP (Thr668)
Concentration:	1 mg/ml
Gene Name:	APP A4 AD1
Protein Name:	Amyloid beta A4 protein, Amyloid- β , A β
Other Name:	Amyloid beta A4 protein (ABPP) (APPI) (APP) (Alzheimer disease amyloid protein) (Cerebral vascular amyloid peptide) (CVAP) (PreA4) (Protease nexin-II) (PN-II) [Cleaved into: N-APP; Soluble APP-alpha (S-APP-alpha); Soluble APP-beta (S-APP-beta); C99; Beta-P05067
Human Swiss Prot No.:	P05067
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Rat APP (phospho-Thr668)
Storage and Stability:	-20°C/1 year