

Antibodies Product Datasheet

Product Name:	CARD11 (Phospho-Ser652) Antibody
Catalog Number:	A9545
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human CARD11 (Ser652)
Concentration:	1 mg/ml
Gene Name:	CARD11 CARMA1
Protein Name:	CARD11 (Ser652)
Other Name:	Caspase recruitment domain-containing protein 11 (CARD-containing MAGUK protein 1) (Carma 1)
Human Swiss Prot No.:	Q9BXL7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Rat CARD11 (phospho-Ser652)
Storage and Stability:	-20°C/1 year