

Antibodies Product Datasheet

Product Name:	MLKL (Phospho-S358) Antibody
Catalog Number:	A9555
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), WB
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human MLKL (phospho S358)
Concentration:	1 mg/ml
Gene Name:	MLKL
Protein Name:	MLKL (phospho S358)
Other Name:	Mixed lineage kinase domain-like protein
Human Swiss Prot No.:	Q8NB16
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous phospho levels of MLKL (phospho S358) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year