

Antibodies Product Datasheet

Product Name:	ARHGAP42 (Phospho-Tyr376) Antibody
Catalog Number:	A9557
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human ARHGAP42 (Tyr376)
Concentration:	1 mg/ml
Gene Name:	ARHGAP42
Protein Name:	ARHGAP42 (Tyr376)
Other Name:	Rho GTPase-activating protein 42 (Rho GTPase-activating protein 10-like) (Rho-type GTPase-activating protein 42)
Human Swiss Prot No.:	A6NI28
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse ARHGAP42 (phospho-Tyr376)
Storage and Stability:	-20°C/1 year