

Antibodies Product Datasheet

Product Name:	Atg14 (Phospho-Ser29) Antibody
Catalog Number:	A9559
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human Atg14 (Ser29)
Concentration:	1 mg/ml
Gene Name:	ATG14 KIAA0831
Protein Name:	Atg14 (Ser29)
Other Name:	Beclin 1-associated autophagy-related key regulator (Barkor) (Autophagy-related protein 14-like protein) (Atg14L)
Human Swiss Prot No.:	Q6ZNE5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Atg14 (phospho-Ser29)
Storage and Stability:	-20°C/1 year