

## Antibodies Product Datasheet

|                               |                                                                                                                                                                                                                          |
|-------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | BACH1/BRIP1 (Phospho-Thr1133) Antibody                                                                                                                                                                                   |
| <b>Catalog Number:</b>        | A9560                                                                                                                                                                                                                    |
| <b>Reactivity:</b>            | Human, Rat, Mouse                                                                                                                                                                                                        |
| <b>Applications:</b>          | WB                                                                                                                                                                                                                       |
| <b>Dilutions:</b>             | WB 1:1000-2000                                                                                                                                                                                                           |
| <b>Modification:</b>          | Phospho Specific                                                                                                                                                                                                         |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                                                                                                                        |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.                                                                                                                |
| <b>Immunogen:</b>             | Synthesized phospho peptide around human BACH1 and BRIP1 (Thr1133)                                                                                                                                                       |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                                                                                                                  |
| <b>Gene Name:</b>             | BRIP1 BACH1 FANCI                                                                                                                                                                                                        |
| <b>Protein Name:</b>          | BACH1/BRIP1 (Thr1133)                                                                                                                                                                                                    |
| <b>Other Name:</b>            | Fanconi anemia group J protein (Protein FANCI) (EC 3.6.4.13) (ATP-dependent RNA helicase BRIP1) (BRCA1-associated C-terminal helicase 1) (BRCA1-interacting protein C-terminal helicase 1) (BRCA1-interacting protein 1) |
| <b>Human Swiss Prot No.:</b>  | Q9BX63                                                                                                                                                                                                                   |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                                  |
| <b>Specificity:</b>           | This antibody detects endogenous levels of Human BACH1 or BRIP1 (phospho-Thr1133)                                                                                                                                        |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                                                                                                             |