

Antibodies Product Datasheet

Product Name:	BAP1 (Phospho-Ser592) Antibody
Catalog Number:	A9561
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human BAP1 (Ser592)
Concentration:	1 mg/ml
Gene Name:	BAP1 KIAA0272 hucep-6
Protein Name:	BAP1 (Ser592)
Other Name:	Ubiquitin carboxyl-terminal hydrolase BAP1 (EC 3.4.19.12) (BRCA1-associated protein 1) (Cerebral protein 6)
Human Swiss Prot No.:	Q92560
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human BAP1 (phospho-Ser592)
Storage and Stability:	-20°C/1 year