

## Antibodies Product Datasheet

|                               |                                                                                                             |
|-------------------------------|-------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | BCL9L (Phospho-Ser915) Antibody                                                                             |
| <b>Catalog Number:</b>        | A9562                                                                                                       |
| <b>Reactivity:</b>            | Human, Rat, Mouse                                                                                           |
| <b>Applications:</b>          | WB                                                                                                          |
| <b>Dilutions:</b>             | WB 1:1000-2000                                                                                              |
| <b>Modification:</b>          | Phospho Specific                                                                                            |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                           |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.   |
| <b>Immunogen:</b>             | Synthesized phospho peptide around human BCL9L (Ser915)                                                     |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                     |
| <b>Gene Name:</b>             | BCL9L DLNB11                                                                                                |
| <b>Protein Name:</b>          | BCL9L (Ser915)                                                                                              |
| <b>Other Name:</b>            | B-cell CLL/lymphoma 9-like protein (B-cell lymphoma 9-like protein)<br>(BCL9-like protein) (Protein BCL9-2) |
| <b>Human Swiss Prot No.:</b>  | Q86UU0                                                                                                      |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                     |
| <b>Specificity:</b>           | This antibody detects endogenous levels of Human BCL9L (phospho-Ser915)                                     |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                |