

Antibodies Product Datasheet

Product Name:	BRSK2 (Phospho-Thr260) Antibody
Catalog Number:	A9563
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human BRSK2 (Phospho Thr260)
Concentration:	1 mg/ml
Gene Name:	BRSK2 C11orf7 PEN11B SADA STK29 HUSSY-12
Protein Name:	BRSK2 (Phospho Thr260)
Other Name:	Serine/threonine-protein kinase BRSK2 (EC 2.7.11.1) (Brain-selective kinase 2) (EC 2.7.11.26) (Brain-specific serine/threonine-protein kinase 2) (BR serine/threonine-protein kinase 2) (Serine/threonine-protein kinase 29) (Serine/threonine-protein kinase S
Human Swiss Prot No.:	Q8IWK3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat BRSK2 (Phospho Thr260)
Storage and Stability:	-20°C/1 year