

Antibodies Product Datasheet

Product Name:	CTDSPL2 (Phospho-Ser104) Antibody
Catalog Number:	A9565
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human CTDSPL2 (Ser104)
Concentration:	1 mg/ml
Gene Name:	CTDSPL2 HSPC058 HSPC129
Protein Name:	CTDSPL2 (Ser104)
Other Name:	CTD small phosphatase-like protein 2 (CTDSP-like 2) (EC 3.1.3.-)
Human Swiss Prot No.:	Q05D32
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human CTDSPL2 (phospho-Ser104)
Storage and Stability:	-20°C/1 year