

Antibodies Product Datasheet

Product Name:	CYLD (Phospho-Ser418) Antibody
Catalog Number:	A9566
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human CYLD (Ser418)
Concentration:	1 mg/ml
Gene Name:	CYLD CYLD1 KIAA0849 HSPC057
Protein Name:	CYLD (Ser418)
Other Name:	Ubiquitin carboxyl-terminal hydrolase CYLD (EC 3.4.19.12) (Deubiquitinating enzyme CYLD) (Ubiquitin thioesterase CYLD) (Ubiquitin-specific-processing protease CYLD)
Human Swiss Prot No.:	Q9NQC7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human CYLD (phospho-Ser418)
Storage and Stability:	-20°C/1 year