

Antibodies Product Datasheet

Product Name:	Desmoplakin (Phospho-Ser165/166) Antibody
Catalog Number:	A9567
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human Desmoplakin (Ser165 and 166)
Concentration:	1 mg/ml
Gene Name:	DSP
Protein Name:	Desmoplakin (Ser165/166)
Other Name:	Desmoplakin (DP) (250/210 kDa paraneoplastic pemphigus antigen)
Human Swiss Prot No.:	P15924
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Rat Desmoplakin (phospho-Ser165 or 166)
Storage and Stability:	-20°C/1 year